



BUSTER Design Collar Bones

Make the sending home of the animal more fun with the new BUSTER Design Collar with printed bones. Watch the surprise and happiness of your clients when the dog wears printed collar in eye-catching colours.


- Quick fastening system
- Colour-coded according to size
- Full body coverage
- Multiple adjustment points for proper fit




Turquoise




Purple




Orange




Green



Yellow



Red



Blue



BUSTER Design Collar Bones

The BUSTER Design Collars are printed collars in eye-catching colours.

The BUSTER Design Collar Bones is transparent and has the quick and simple clic closing system, which is well- known from the BUSTER Clic Collars. The collars are made of 100% polypropylene.

Unique features and benefits

Safety and perfect fit

The collar is provided with multiple adjustment points for proper and safe fit.

Easy application

The collar has a quick and simple clic closing system, which makes the application of the collar as easy as buttoning your shirt.

Colour-coded according to size

The collar is a printed collar with bones in eye-catching colours; one colour for each size of collar for easy recognition.

Full body coverage

The collar offers a full body coverage and protection preventing the animal from reaching certain areas of the body.

Recommended use for BUSTER Design Collar Bones

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and other skin ailments
- Skin diseases

Instructions for use


Scan QR code to see the video



shorturl.at/IDKU1



Cat. No	Size	Pk. Size	Colour
273731	7.5 cm	10/pk	
273732	10.0 cm	10/pk	
273733	12.5 cm	10/pk	
273734	15.0 cm	10/pk	
273735	20.0 cm	10/pk	
273736	25.0 cm	10/pk	
273737	30.0 cm	10/pk	



Size	Length in cm / Inch	Opening in cm / Inch	Neck-hole in cm / Inch	Suitable for	Lbs
7.5 cm	7.5 cm / 2.95"	13 - 15 cm / 5.1" - 5.9"	7 - 8 cm / 2.7" - 3.1"	Pomeranian, Chihuahua	1-6
10 cm	10 cm / 3.93"	15 - 18 cm / 5.9" - 7.0"	7 - 8 cm / 2.7" - 3.1"	Yorkshire terrier, Bichon Frise, Maltese	6-10
12.5 cm	12.5 cm / 4.92"	19 - 21cm / 7.4" - 8.2"	9 - 10.5 cm / 3.5" - 4.1"	Danish/Swedish farmdog, Papillon, Cairn Terrier, Border Terrier	10-20
15 cm	15 cm / 5.90"	23 - 27 cm / 9.0" - 10.60"	10 - 12 cm / 3.9" - 4.7"	Dutch Kooikerhondje, Scottish Terrier, Beagle, Irish Terrier	20-30
20 cm	20 cm / 8.87"	20 - 22 cm / 7.8" - 8.6"	12 - 14 cm / 4.7" - 5.5"	Australian Shepherd, English Springer Spaniel, Staffordshire Bull Terrier	32-50
25 cm	25 cm / 9.84"	34 - 42 cm / 13.3" - 16.5"	14 - 16 cm / 5.5" - 6.2"	Labrador Retriever, Siberian Husky, Irish Setter	50-60
30 cm	30 cm / 11.81"	41 - 48 cm / 16.1" - 18.8"	15 - 18 cm / 5.9" - 7.0"	Golden Retriever, Boxer, Collie, German Shepherd	50-70

