

# THANK YOU FOR YOUR PURCHASE!

“When a feline patient bit me early in my veterinary career, my injury was potentially life-threatening. I developed **ArmOR Hand Protective Gloves** to protect animal caregivers from painful injuries that can be both expensive and traumatic.



With your purchase of **ArmOR Hand**, you are supporting the *One Health Mission*, which emphasizes the undeniable connection of humans, animals, and our environment.

Thank you for the work that you do,

*Laura Catena*  
LAURA CATENA, DVM

THE  
**ARMORHAND**  
PROTECTOR

[ARMORHANDGLOVE.COM](http://ARMORHANDGLOVE.COM)



THE  
**ARMORHAND**  
PROTECTOR

**MORE FEELS. LESS FORCE.**

Comfortable, form-fitting, & durable  
A must-have for animal handlers

## DETERMINE YOUR SIZE USING THE GUIDELINES BELOW:

2X-Large 4.0" - 4.25"  
(102-108mm)  
X-Large 3.75" - 4.0"  
(95-102mm)  
Large 3.5" to 3.75"  
(90-95mm)  
Medium 3.25" to 3.5"  
(83-90mm)  
Small 3.0" to 3.25"  
(76-83mm)  
X-Small 2.75" to 3.0"  
(70-76mm)

PLEASE USE A CLOTH TAPE MEASURE  
ON THE HAND YOU USE THE MOST

THERE ARE TWO MEASUREMENTS TO  
DETERMINE YOUR SIZE:

1. LENGTH OF MIDDLE FINGER
  - A. X-Small 2.75" to 3.0" (70-76mm)
  - B. Small 3.0" to 3.25" (76-83mm)
  - C. Medium 3.25" to 3.5" (83-90mm)
  - D. Large 3.5" to 3.75" (90-95mm)
  - E. X-Large 3.75" - 4.0" (95-102mm)
  - F. 2X-Large 4.0" - 4.25" (102-108mm)
2. CIRCUMFERENCE OF PALM (EXCLUDE THUMB)
  - A. X-Small = 7.5" to 8.0" (191-203mm)
  - B. Small = 8.0" to 8.5" (203-216mm)
  - C. MEDIUM = 8.5" to 9.0" (216-228mm)
  - D. LARGE = 9.0" to 9.5" (228-240mm)
  - E. X-LARGE = 9.5" to 10.0" (240-254mm)
  - F. 2X-LARGE = 10.0" to 10.5" (254-267mm)

IF YOU FALL BETWEEN SIZES — OPT FOR THE LARGER SIZE

\*\*\* DIAGRAM NOT TO SCALE \*\*\*

